



R-4608 Series

aluminum 4-wheeled rollators w/8" wheels

- 4 wheeled, folding, lightweight aluminum rollator (rolling walker) with powder coat finish in optional colors
- Adjustable height, Hand-grips with locking loop brakes on 8" casters
- Padded large seat and curved backrest for added comfort
- Easy-remove storage basket with handle

	seat pad Width (in.)	seat pad depth (in.)	seat width (in.)	seat height (in.)	backrest height (in.)	overall height (in.)	overall open width (in.)	weight (lb.)	weight capacity (lb.)
R-4608	14	15	18	23	11	35~40	24	15	250

Model No.	Description
R-4608	4-wheel rollator with loop (grip) brakes, seat and basket

Color: BL=blue, BD=burgundy

** Price & specifications are subject to change without notice



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335